

Welcome to the 12th Annual MAVES

3 May 2020



Thank YOU, 2020 vendors, for your support!!

During the breaks, please visit our supporting vendors in their breakout room through the CHAT function and be entered into a raffle for Amazon gift cards!

BOEHRINGER



PLATINUM SPONSORS

Caring for Lives through Innovation, Quality and Service

LivaNova

Health innovation that matters

Spectrum Medical

TERUMO



Invosurg


SILVER SPONSORS



Abbott

Extracorporeal Membrane Oxygenation: Just the basics

Sarah Guthrie, PA-C



"Copyrighted material used with permission of the author, University of Iowa Hospitals & Clinics, uihc.org."



Objectives

- ❑ Understand the basics of ECMO and its application
- ❑ Identify key components to the ECMO circuit
- ❑ Identify the role of Clinical Specialist in the care of ECMO patient
- ❑ Identify methods for anticoagulation
- ❑ Recognize additional mechanical interventions seen with ECMO

ECMO basics

- Used for management of life threatening pulmonary and/or cardiac failure when other treatment is not working
- Essentially a modification of the cardiopulmonary bypass circuit used in cardiac surgery
- Temporary support – allow for time for treatment and recovery of organs
- Used to deliver oxygen and remove carbon dioxide
- Two types:
 - Respiratory Failure requires VV (veno-venous) ECMO configuration
 - Cardiac Failure requires VA (veno-arterial) ECMO configuration

What is Extracorporeal Membrane Oxygenation?

- As defined by the Extracorporeal Life Support Organization (ELSO):

The use of mechanical devices to support heart and/or lung function in severe heart or lung failure, unresponsive to optimal conventional care. (Brogan, 2017)



Goals of ECMO

- Bridge to decision
 - Recovery
 - Durable ventricular Assist Device
 - Transplant
 - Withdrawal

Indications for ECMO

Respiratory Failure (V-V)

- Refractory Hypoxemia
- Refractory Hypercapnia
- Ventilator Induced Lung Injury
- Acute Respiratory Distress Syndrome
- H₁N₁

Cardiovascular Compromise (V-A)

- Acute myocardial infarction
- Post cardiectomy
- Acute Myocarditis
- Acute Pulmonary Embolism
- Bailout post cardiac intervention
- Preoperative Support
- Acutely decompensated chronic cardiomyopathy

Exclusion criteria of ECMO

Absolute Contraindication

- Pre-existing condition incompatible with recover
 - Severe neurological injury
 - End stage malignancy

Relative Contraindication

- Very poor prognosis from primary condition
- Uncontrollable bleeding
 - Unable to tolerate anticoagulation
- Aortic insufficiency
- Aortic Dissection




Types of ECMO support access

- Percutaneous veno-venous (V-V) ECMO
- Central veno-atrial (V-A) ECMO
- Percutaneous veno-atrial (V-A) ECMO
- Percutaneous veno-venous-atrial (V-V-A) ECMO
- Percutaneous veno-atrial-venous (V-A-V) ECMO

Percutaneous V-V ECMO

Traditional cannulation placement

- Drainage cannula- right femoral vein, tip in IVC
- Right internal jugular, tip in SVC
- Can also cannulate via bilateral femoral veins



Percutaneous V-V ECMO

Dual Cannula

- drain from distal port in IVC
- return flow into RA




Figure 3: single cannula V-V ECMO (Napp et al, 2016)

Central V-A ECMO

- Drainage cannula in right atrium
- Return cannula in aorta
- Requires sternum to remain open
- Difficult to transport





Figure 1: Central V-A ECMO ([Pavlushkov et al., 2017](#))

Percutaneous V-A ECMO


- Drainage cannula in femoral vein terminating by RA
- Return cannula in femoral artery, terminating in Iliac/aorta
- Need for flow into the SFA of cannulated limb
 - Antegrade perfusion canula- Return flow cannulas to prevent limb ischemia



Percutaneous V-V-A ECMO

Reduces circulatory overload of the heart and reduces filling pressures


- Blood is drained from SVC and IVC
- Further lowers blood volume into the heart/lungs
- Effectively decreases CVP and PA pressures



Percutaneous V-A-V ECMO

Return flow is divided

- Drainage cannula in femoral vein, terminating in IVC by RA
- Return cannulas in femoral artery, terminating in iliac/aorta, and right internal vein, terminating in SVC by the RA



Equipment

- Perfusionist typically maintain the ECMO circuits and their availability
 - cannulas
 - tubing
 - pump
 - membrane oxygenator
 - heat exchanger
 - gas blender
- Each facility maintains their own set up and device that varies in size and transport capability

Basic ECMO Circuit

1. Pump

- RPM
- Volume
- Flow

2. Membrane Oxygenator

- FiO_2
- Sweep

3. Heat Exchanger




Figure 7: Circuit ([emDOCs, 2018](#))



Pumps

Centrifugal

- More efficient designs of pumps result in safer use than previous designs
- Magnetic driven and magnetic suspended pumps available

Roller

- Positive displacement pump generating forward flow as a function of tubing size and pump speed
- Fundamental in the conduct of CPB and ECLS
- Better for lower blood flows



<https://thoracickey.com/extracorporeal-membrane-oxygenation-in-infants-and-children/>



Flow

- V-A ECMO (for cardiac failure): 50-60cc/kg/min. The flow is limited by vascular access, length and size of cannula, and pump properties.
- V-V ECMO (for respiratory failure): 60-80cc/kg/min. CO₂ removal always exceeds O₂ delivery. Besides blood flow, oxygenator/membrane properties and gas sweep determine O₂ and CO₂ levels.

Membrane Oxygenator

Used to add oxygen and remove CO₂

Sweep: gas ventilated through gas exchange



Membrane Oxygenator

ABG goals on ECMO

- pH: 7.35-7.45
- pCO₂: 35-45 mmHg
- pO₂: 80-100mmHg

Titrating FiO₂:

- $\uparrow \text{FiO}_2 = \uparrow \text{pO}_2$
- $\downarrow \text{FiO}_2 = \downarrow \text{pO}_2$

Titrating Sweep:

- $\uparrow \text{sweep} = \uparrow \text{CO}_2 \text{ removal}$
- $\downarrow \text{sweep} = \downarrow \text{CO}_2 \text{ removal}$

Delivered Oxygen

- $\text{DO}_2 = \text{CO} \times \text{CaO}_2$ (mL/min/m²)
- $\text{CaO}_2 = (1.34 \times \text{Hgb} \times \text{SaO}_2) + (0.003 \times \text{PaO}_2)$

Clinical Specialists

- ECMO clinical specialist (CS) is "technical specialist trained to manage the ECMO system and the clinical needs of the patient"
- CS may be perfusionists, nurses, respiratory therapists
 - Perfusionists may be ideal due to their knowledge and training of CPB and ECMO
- Defined CS help streamline and provide consistency for ECMO emergencies
- Resource for the providers, nurses, patient



Anticoagulation Strategies

- Unfractionated Heparin (UNFH) infusion
 - First line medication for ECMO thrombus prevention
 - Measured by aPTT
 - Risk to develop heparin induced thrombocytosis and thrombosis (HITT)
 - Risk to develop heparin resistance

Anticoagulation Strategies

- Bivalirudin (Angiomax) infusion
 - DTI- direct thrombin inhibitor
 - Measured by aPTT
 - Ideal in setting of HIT, heparin resistance, or non-HIT thrombocytopenia
 - Falsely raises INR
 - Renally cleared, can be used in CRRT/HD

Anticoagulation Strategies

- Argatroban infusion
 - DTI
 - Measured by aPTT
 - Second line medication when bivalirudin is not ideal
 - Falsely raises INR
 - Contraindicated in transaminitis/shock liver



Anticoagulation Monitoring

- Activated Clotting Time (ACT)
- Activated Partial Thromboplastin Time (aPTT)
- Anti-factor Xa Assay (Anti-Xa)
- Antithrombin Level (AT)
- Lactic Dehydrogenase (LDH)

ACT: goal 160-180 seconds

Advantages

- Can be done at bedside
- Minutes to obtain
- Only a drop of blood

Disadvantages

- Multiple causes for high ACT
 - Excessive anticoagulation
 - Thrombocytopenia
 - Coagulopathy
 - Combination of all above
- Variability in the ACT even from a single sample

aPTT goal: 60-80 seconds

Advantages

- More accurate than ACT
- Decreased risk of hemorrhagic complications when using UNFH

Disadvantages

- Increased risk of circuit clots when using UNFH
- Unreliable in critical illness due to effects of acute phase reactants
 - Falsely prolonged with elevated C-reactive protein
 - Falsely decreased with elevated Factor VIII

Anti-Xa goal: 0.3-0.5 units/mL

Advantages

- More accurate than aPTT
- Decreased blood product use, hemorrhagic complications, and increased circuit life

Disadvantages

- Send out lab at most hospitals – prolonged turn around time and costs

Antithrombin [deficiency]

- Antithrombin is an Important inhibitor of coagulation
 - need it to stop coagulation
- AT deficiency occurs secondary to chest tube losses (consumption)
- Can be replaced by fresh frozen plasma or by Antithrombin III (AT₃)
 - Counterintuitive– need to restore antithrombin to prevent clots
- Send out lab at most hospitals



LDH and Plasma Free Hemoglobin

- Measures of hemolysis from the ECMO circuit due to:
 - Sheering stress
 - High ECMO flow
 - Cavitation (chugging, chatter)
 - Pressure changes in oxygenator
- Elevated values increase risk of circuit thrombus and embolic events

Lab Draws

- How often? 4 hours? 6 hours? 12 hours?
 - ABG
 - Arterial blood gas – bedside (EPOC, iStat)- quick result, less accurate, more expensive
 - Conventional – lab – accurate (gold standard), slower results
 - CBC
 - PT/INR/PTT
 - Fibrinogen
 - LDH
 - Lactic Acid

To transfuse or not?

| Lab evaluation | Goal |
|----------------|--|
| Platelets | Transfuse to maintain $> 100,000$ uL |
| INR | FFP transfusion to maintain $INR < 2$ |
| Fibrinogen | Cryoprecipitate to maintain fibrinogen > 100 mg/dL |
| Hematocrit | PRBCs to maintain hematocrit $> 25\%$ |
| Antithrombin | FFP or AT ₃ to maintain Antithrombin $> 50\%$ |

Pulsatile vs Non-pulsatile management

Pulsatile


- SBP goal 100-120mmHg
- LV contracting reduces risk of developing a thrombus
- Potential for North-South Syndrome

Non-pulsatile

- MAP goal 60-80 mmHg
- May need surgical LV vent or Impella to decompress and prevent thrombus formation

North-South Syndrome

- Aka: Harlequin Syndrome, “Blue Head, Red Legs”
- Recovering heart prevents the retrograde flow from the ECMO circuit to perfuse the upper half of the body
- Consider V-A-V ECMO if unable to decannulate



Left Ventricle Vent

Reasons

- Findings of LV dysfunction with increased LV end-diastolic and systolic volumes increase myocardial wall stress
- Ex: impaired LV unloading, LV stasis, and pulmonary edema
- Is it needed?

Methods


- Surgical vent – cannula placed in left atrium or ventricle or pulmonary artery
 - Axillary artery through the aortic valve to ventricle
- Intra-aortic balloon pump-decompresses by lowering afterload reduction
- Impella – effectively drains LV

CRRT and ECMO

A: First-line connection of CVVH to the ECMO circuit. Both lines connected between pump and oxygenator (segment B).

B: Second-line connection: Access line connected before the pump (segment A).

C: Both lines before the pump (segment A). Optional clamp adjusted on return line if pressure is below the low-pressure alarm. CVVH, continuous veno-venous hemofiltration.



Option 1: Pressures in CVVH access line < +200 mmHg and in return line < +350 mmHg

Advantages:

- No risk of air embolism
- Return line connected on positive pressure
- No oxygenation membrane shunt
- No recirculation

Drawbacks:

- Access line connected on positive pressure

Option 2a: Used for high positive pressures in CVVH access line: > 200 mmHg

Advantages:

- Access line connected on negative pressure
- Return line connected on positive pressure
- No oxygenation membrane shunt
- No recirculation

Drawbacks:

- hypothetical risk of air embolism

Option 2b: Used for high positive pressures in CVVH return line: > 350 mmHg

Advantages:

- Access line connected on negative pressure
- No oxygenation membrane shunt
- No recirculation

Drawbacks:

- Theoretical risk of air embolism
- Return line connected on negative pressure
 - But Clamp can mitigate negative pressure on return line in case of pressure below low pressure alarm

S.A: Segment A: ECMO segment with negative pressure

S.B: Segment B: ECMO segment with positive pressure

Patient Management - Ventilator

- Central V-A ECMO
 - Rest settings (low VT, PEEP support)
 - Goal is to maintain lung volumes
- Percutaneous ECMO
 - Variable, modes/settings are dependent on clinical situation and the presence of North-South Syndrome
- No evidence basis for ideal strategy, from anecdotal experience



Patient Management- Vasopressors

- Epinephrine, vasopressin, norepinephrine, phenylephrine
- Extremity ischemia/compartiment syndrome



Patient Management - Inotropes

- Milrinone, epinephrine, dopamine, dobutamine
- No evidence basis for idea strategy, from anecdotal experience

References

- Annich G.M., Lynch W.R, MacLaren, G., Wilson J.M., Bartlett, R.H. (2012). *ECMO Extracorporeal Cardiopulmonary Support in Critical Care* (4th ed). Ann Arbor, MI: ELSO
- Brogan T.V., Lequier L. Rovertto L., Graem M., Peek G. (2017). *Extracorporeal Life Support: The ELSO Red Book* (5th ed). Ann Arbor, MI: ELSO
- Camboni D, Schmid C. To vent or not on veno-arterial extracorporeal membrane oxygenation, does it improve myocardial recovery and outcome?. *J Thorac Dis.* 2017;9(12):4915–4918. doi:10.21037/jtd.2017.11.98
- de Tymowski, Christian; Augustin, Pascal; Houissa, Hamda; Allou, Nicolas; Montravers, Philippe; Delzongle, Alienor; Pellenc, Quentin; Desmard, Mathieu *ASAIO Journal*63(1):48-52, January/February 2017. doi: 10.1097/MAT.0000000000000441
- Pavlushkov E, Berman M, Valchanov K. Cannulation techniques for extracorporeal life support. *Ann Transl Med* 2017; 5(4): 70. doi: 10.21037/atm.2016.11.47
- S Strunina, P Ostadal. Left ventricle unloading during veno-arterial extracorporeal membrane oxygenation. *Curr Res Cardiol* 2016;3(1):5-8.
- emDOCs [Image]. 2018. Retrieved from <http://www.emdocs.net/extracorporeal-membrane-oxygenation-ecmo-in-the-ed/>